

Methoxyflurane (Penthrox®)

	Ensure the Activated Carbon (AC) Chamber is inserted into the dilutor hole on the top of the inhaler.
	2. Remove the cap of the bottle by hand. Alternatively, use the base of the inhaler to loosen the cap with a ½ turn. Separate the inhaler from the bottle and remove the cap by hand.
	3. Tilt the inhaler to a 45° angle and pour the total contents of one bottle into the base of the inhaler whilst rotating.
	4. Place wrist loop over patient's wrist. Patient inhales through the mouthpiece of the inhaler to obtain analgesia. First few breaths should be gentle and then breathe normally through Inhaler.
5	5. Patient exhales into the inhaler (to protect you). The exhaled vapour passes through the AC Chamber to adsorb any exhaled methoxyflurane.
	If stronger analgesia is required, patient can cover dilutor hole on the AC Chamber with finger during use.
	7. If further pain relief is required, after the first bottle has been used use a second bottle if available. Alternatively use a second bottle from a new combination pack. Use in the same way as the first bottle in Step 2 and 3. No need to remove the AC Chamber. Put used bottle into the plastic bag provided. Patient should be instructed to inhale intermittently to achieve adequate analgesia. Continuous inhalation will reduce duration of use. Minimum dose to achieve analgesia should be administered.
	8. Replace cap onto bottle. Place used inhaler and used bottle in sealed plastic bag and dispose of responsibly in a designated bin or suitable waste container.